INSTRUCTIONS FOR USE

CyberCement Resin

Permanent adhesive resin cement - dual cured

INGREDIENTS:

Base & Catalyst:

- 4-Methacryloxyethyltrimellitic anhydride.
- Urethan dimethacrylate.
- Multifunctional monomers.
- Polymerization activators.
- Photoinitiators.
- Glass filler.

INDICATIONS:

Permanent cementation of crowns and bridges, inlays and onlays, posts and cores, ceramic crowns and Maryland bridges, veneering of alloys.

PROPERTIES:

CyberCement Resin is a multi-purpose, dual cured, radiopaque, permanent cement. It can be used with base alloys, porcelain, dentine, enamel and with amalgams. Light curing of CyberCement Resin gives a better strength at the margins and superior stability.

SETTING TIME:

Working time at room temperature approximately 2 minutes. Oral set time is approximately 3,5 minutes. Light cure for 20 seconds.

CAUTION:

CyberCement Resin contains polymerisable monomers which may cause allergic reactions with susceptible persons. Wash thoroughly with soap and water after contact. If skin sensitivity occurs or if known allergy to methacrylate resin exists, discontinue use.

STORAGE

- Do not store above 24 °C (75 °F).
- Keep away from direct sunlight.
- For optimum freshness, keep refrigerated.
- Do not freeze.
- If refrigerated, allow the syringe to reach room temperature.
- Do not use after expiration date.See expiry date on the label.
- Shelf life is 2 years after the date of manufacturing.

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CB Healthcare Consulting GmbH